




NVoy Polymer - General

What is it?




- ④ Specially chosen carbohydrate backbone with selected hydrophobic side chains derivitised at specific positions
- ④ Linear in solution with side chains aligning to form Hydrophobic face
- ④ $M_w = 5\text{kDa}$ ($R_{\text{hyd}} \sim 18\text{kDa}$), Does not access binding sites
- ④ Association very dynamic, NVoy Polymer is very easy to completely remove from sample
 $k_d = 28 \mu\text{M}$ (GFP) & $12 \mu\text{M}$ (Hexokinase) – See Removal Protocol for more information
- ④ NVoy is Uncharged, UV transparent and pH Stable (pH 2 to 11)
- ④ Does not affect existing purification methodologies




NVoy Polymer - General

- NVoy works by associating with any surface exposed hydrophobicity, masking it and presenting the hydrophobic backbone to the solvent.
- This prevents hydrophobic interaction between the target protein and 1) other target proteins preventing aggregation and promoting a heterogeneous sample, 2) process surfaces (i.e. purification resin, filtration membranes) increasing yields, 3) contaminant proteins improving purity of sample, reducing aggregation and increasing yields
- Hydrophobic interaction can be cited as the cause for protein aggregation in ~75% of cases (Key Note presentation at PEGS 2007 meeting)






Aggregation Suppression





Replacing Detergents

- Preparation of a cytosolic enzyme involved in DNA synthesis
- Enzyme required for screening and structural biology but normal production in detergents was not suitable for crystallography
- Protein processed in and purified with
- NVoy polymer in buffers gave increase in both the purity and yield of the protein.
- Specific activity of protein was similar to literature values for native enzyme







Replacing Detergents

- Can be used to completely replace detergents in membrane protein preparation
- Optimal NVoy concentration for extraction is 15 mg/ml although up to 25 mg/ml has been used.

Histidine Kinase Receptor




E. coli IMP tagged with GFP





Replacing Detergents




- On column detergent replacement* performed using IMP-GFP and Histidine Kinase extracted in FC12
- Analysed by NMR – All FC12 replaced by NVoy
- Confirmed that NVoy analogous to detergent mode of action





Cell Free Expression

- Family A and B GPCR's expressed CF Conventionally and with NVoy
- Found that NVoy works as in Mode B
- Concluded that NVoy does not interfere with CF expression
- NVoy enables soluble expression of active GPCR's (>10 nM ligand affinity)





Increase Purity & Stability

- Large membrane associated protein (>160kDa) that requires co-factors for activity


Using Detergents;

results in unstable enzyme preparations containing large number of contaminants

- At 4°C loss of activity within 24h
- At -80°C after one freeze/thaw cycle

Using NVoy Polymer;

- Cleaner protein preparations obtained
- Enzyme activity stable at 4°C for over 4 weeks
- Enzyme activity stable at -80°C for 4 freeze/thaw cycles



- Enzyme processed from *E. coli* lysates
- Cells cultured at 37°C, induced with IPTG and cultured for further 72h.
- Cells harvested then lysed with buffer containing co-factors and NVoy polymer / detergents Lysate clarified by centrifugation and purified on ADP Sepharose using NVoy polymer



Remove Endotoxin

- A method that can be introduced into any existing purification process
- Comparison of three methods
 - Pierce Detoxi-gel
 - Non-ionic detergent (Triton X114) according to Reichelt *et al.*, 2006
 - NV polymer

Endotoxin removal - His tagged protein was bound to Talon resin

MBP fusion protein was bound to MBP resin

3 washing steps were carried out with the protein held on the affinity columns

1. 10CV of buffer containing 2% NV10
2. 10CV of 0.2% NV10
3. 5 CV of buffer at 4°C

Method	Protein Yields (% Control)	Endotoxin EU	Activity (% Control)
Control MBP	100	120,000	100
Pierce	32	3000	55
Triton X114	45	80	32
NVoy Polymer	78	None Detected	72
Control His	100	130,000	100
Pierce	29	1200	37
Triton X114	39	100	25
NVoy Polymer	66	<1	77



Improve Concentration

BSA Solutions concentrated tenfold, as per manufacturers instructions, using 5,000 MWCO VivaSpin 2 spin column. Addition of NVoy Polymer dramatically increased the recovered yield by preventing the non-specific binding to the HydroSart filtration membrane.

Starting solution (1 ml)	Recovered yield (%)
10mg/ml BSA	46%
10mg/ml BSA + 10mg/ml NVoy	60%
10mg/ml BSA + 40mg/ml NVoy	85%
10mg/ml BSA + 100mg/ml NVoy	90%

This type of experiment has been performed by numerous people using different proteins with excellent results in all cases to date.



Improve Concentration

- Sample of MrB – a mouse murine morphogen binding protein and a type I membrane protein which is bound to the extracellular membrane.
- MrB is Important because it regulates a pathway involved in specifying cell differentiation during embryogenesis; dysfunction of MrB is associated with cancers and congenital defects.
- Concentration prior to crystallisation – 60% loss!

Sample, MrB +	[Final] (mg/ml)	Final Volume (μ l)	Amount Recovered (μ g)	Yield (%)
0 mg/mg NVoy	3.5	20	70	44.2
0.05 mg/mg NVoy	4.4	25	110	69.4
0.25 mg/mg NVoy	5.2	26	135.2	85.4
1 mg/mg NVoy	9.7	16	155.2	98.0

