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SCIENCE & TECHNOLOGY

NO-PREP MASS SPEC

New atmospheric pressure ionization method is suitable for security and defense applications

CELIA M. HENRY, CAEN, WASHINGTON

RONNIE A. COOK, the mass spectrometry product manager at JEOL USA in Frisco, Mass., has had one of his last year's projects in mind: a mass spectrometer that can analyze a sample in real time. Just by sticking samples in front of a mass spectrometer, he analyzed paper money, chocolate, hot peppers, and even a needle-laced with a trace of nitrogllycerin.

He is able to do this no-sample preparation method that he and James A. Laramore of GAI Corp., Abingdon, Md., invented. They call the method DART, for direct analysis in real time. DART, which is currently being developed by JEOL, uses a spray of highly charged ions to ionize the sample, which is then analyzed by the mass spectrometer. The method is suitable for security and defense applications.

DART works by applying an electrical potential to a gas such as nitrogen or helium to form a plasma of excited atoms and molecules that then interact with the sample and the atmosphere. Different ionization mechanisms are possible, and operating conditions can be manipulated to favor one over the others. For example, proton transfer is the dominant mechanism of positive ionization.

This type of ionization occurs when molecules in the atmosphere react with water in the atmosphere to produce oxidized water clusters that can penetrate the sample molecule, forming positively charged ions. Under different conditions, electrons

with atmospheric water to form negative-ion clusters of oxygen and water that in turn react with analytes to form negatively charged ions.

In the negative-ionization mode, in turn, and the ions are not protonated because, in DART, plasma formation from the carrier gas is isolated from the air. These ions can interfere with the detection of nitrogen-based explosives and reduce the sensitivity of anion detection.

Cook described the technique publicly for the first time in late January at the 17th National Conference on Mass Spectrometry, John D. (Jack) Henion, a mass spectrometry who serves as chairman, president, and chief executive officer of Johnson BioSciences, Ithaca, N.Y., saw that presentation.

"IT APPEARED as a really simple and easy to use and provided easily interpretable data," Henion told CAEN. "Unanswered questions include whether it will handle larger molecules like peptides and proteins, and can be used for LC/MS applications." However, for the kinds of applications Cook described, it appears like a simple and exciting ionization technique near the commercialization and implementation and new applications.

"JEOL is so excited about the prospects for DART that it created a spin-off company in October 2004 called IonSense to explore the possibility of partnering with manufacturers of other detection technologies. Because the mandate of JEOL's U.S. office doesn't include product development, IonSense will carry out any further development of DART."

The DART intellectual property has been licensed to IonSense, but for now DART will be sold exclusively to JEOL for interfacing with ACQUITY, says Michael Reuter, vice president and transmission electron microscope and scanning probe microscopy product manager at JEOL. IonSense will commercialize the technology for a variety of applications, not just mass spectrometry.

For more, IonSense is a start-up company based at JEOL's headquarters. Members of the DART team include Celia Kerker, and J. Douglas Meibohm, analytical instrument director at JEOL. At this point, none of them is leaving his current position. ■

ANALYZE THIS DART co-inventor Cook shows how easy it is to prepare a sample for analysis: Just stick it in front of the spectrometer.

can be formed if the carrier gas can form ionizable species with high enough internal energy. For example, helium reacts

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


C


D


Joanne Y. Yew, Robert B. Cody, and Edward A. Kravitz, Cuticular hydrocarbon analysis of an awake behaving fly using direct analysis in real-time time-of-flight mass spectrometry, *PNAS*, **105**, 7135-7140 (2008)

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






Waters

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
DART-TOF MS
medium resolution (1~ 6000 fwhm)


DART-Orbitrap MS
high resolution (10k-100k fwhm)

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- Robert B. Cody, Direct Analysis Under Ambient Conditions without Sample Preparation: "Instant Gratification" for Small Molecule Analysis, Conference on Small Molecule Science **CoSMoS 2005** Presentations. August 8-11, 2005
- Harris, G. A.; Fernandez, F. M. Simulations and Experimental Investigation of Atmospheric Transport in an Ambient Metastable-Induced Chemical Ionization Source. *Anal. Chem.*, **81**, 322-329, 2009
- Jana Hajslova, Lukas Vaclavik, Lipodomic fingerprinting: a novel approach to detect fraudulent fats containing high levels of contaminants (3-MCPD esters and others), Detection Technologies for Intentionally Added Adulterants/Contaminants, 11th Annual Fera/JIFSAN Symposium, York – UK, 16-18 June 2010
- Joseph Tice, Douglas Simmons, Michael Festa, James A. Hill, Brian D. Musselman, Development of a "Greener" Ionization source for Ambient Desorption Ionization with Nitrogen Gas, *ASMS 2009*,
- Jana Hajslova, The Use of DART MS for Characterizing Food, 10th joint CSL / JIFSAN symposium on Food Safety and Nutrition, Methods and Systems for Tracking, Tracing, and Verifying Foods, MD, USA, May 13-15 2009,
- Matthew Sharman, Richard J Fussell & Danny Chan, Ambient ionization coupled to hybrid MS: A new technique for food quality / safety control?, 124th AOAC Annual Meeting & Exposition, Florida **USA**, September 26-29, 2010
- Justin Gordon; Elizabeth Crawford; Jing-Tao Wu; Brian D. Musselman; Ming-xiang Liao; Bei-Ching Chuang; Cindy Xia; David Ho; Lily Li; Shaoxia Yu, Quantitative Analysis of Dried Blood Spots by DART (Direct Analysis in Real Time) /MS/MS without Sample Preparation, 58th ASMS Conference on Mass Spectrometry (ASMS 2010), Utah, USA

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- Robert B. Cody, DART and Trace Evidence (a bag of tricks for using DART), National Forensic Science Technology Center Largo, Florida, August 12-13, 2008
- Luke K Ackerman, Gregory O Noonan, Timothy H Begley, 4th International Symposium on Food Packaging, Praha, CzR November 19-21, 2008
- Screening of Forensic Drug Samples using Direct Analysis in Real Time (DART) with Time of Flight Mass Spectrometry :Pamela Reynolds Jeff Leibowitz and Robert Mothershead Chemistry Unit FBI Laboratory
- E.Jagerdeo, J. Leibowitz, J. Clark, L. Schumacher, M. Lebeau, Rapid Analysis with DART™ Source Coupled to TOF; Forensic Applications, Pittcon 2006, Orlando
- Joanne Y. Yew, Robert B. Cody, and Edward A. Kravitz, Cuticular hydrocarbon analysis of an awake behaving fly using direct analysis in real-time time-of-flight mass spectrometry, *PNAS*, **105**, 7135-7140 (2008)
- 草井明彦、リアルタイム直接質量分析法の原理と応用、ぶんせき、124 (2007)