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

The screenshot displays the 'Application Editor' window for 'LOreal Formulation F1.app'. The interface includes a menu bar (File, Edit, Application, Tools, View, Help), a toolbar, and a 'Tasks' panel on the left with icons for various operations like Microwave, Transport task, Decap, Cap / Crimp, Control Cryostat, Transfer Solid, Transfer Liquid, Wait, Set Timer, Switch Contact, Vortex, Heat / Cool, Reflux, Set Vacuum, Set Drawer Valve, Control pH, Set Feed Rate, and Homogenize. The main workspace shows a task list with columns for Task, Name, Parameter, and Description. The selected task is '1 Transfer Liquid' with parameters for liquid transfer from F1MP3 to Formulation Vessel. Below the task list, the 'Source and Destination' configuration panel is open, showing settings for Source Zone (F1MP3), Destination Zone (Formulation Vessel), Equilibration time (0s), Needle Distance (65mm), and Needle configuration (Needle 1 checked).

Task	Name	Parameter	Description
1	Transfer Liquid	Liquid Transfer from Ethanol 1 mL to Waste Syringe 1	Loop Solvent Et
2	Wait	Waiting for 1:00:00 hours	
3	Macro Task	Macro Task	Premix 3 F5
1	Transfer Liquid	Liquid Transfer from F1MP3 to Formulation Vessel	
2	Stir	Agitation ON on zone Formulation Vessel	
3	Set Feed Rate	Start dispense on zone "Feed Pump" from Port C to Port A1 and wait for it to finish	
4	Heat / Cool	Thermostat ON on zone Formulation Vessel	
5	Stir	Agitation ON on zone Formulation Vessel	
6	Transfer Solid	Solid Transfer from Solid F1MP4 to F1MP4	
7	Transfer Solid	Solid Transfer from Solid F1MP5 to F1MP5	
8	Wait	Waiting for 45:00 minutes	
9	Heat / Cool	Thermostat ON on zone Formulation Vessel	
10	Transfer Liquid	Liquid Transfer from Water 1mL to Waste Syringe 1	Loop Solvent W
11	Stir	Agitation ON on zone Formulation Vessel	
12	Transfer Liquid	Liquid Transfer from Phase B to Formulation Vessel	
13	Heat / Cool	Thermostat ON on zone Formulation Vessel	
14	Stir	Agitation ON on zone Formulation Vessel	
15	Wait	Waiting for 30:00 minutes	Addition of F1M
16	Heat / Cool	Thermostat ON on zone Formulation Vessel	
17	Wait	Waiting for 5:00 minutes	

FORMAX[®]


Chemspeed Technologies offers hardware and software for every step in the product development cycle

FORMAX robotic platform

exploded view of the formulation vessel

selection of exchangeable robotic tools

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