



# Newomics M3 Emitters for Shotgun Lipidomics

## Authors

Maoyin Li, Pan Mao,  
Daojing Wang  
Newomics Inc.  
Berkeley, CA 94710

## Goal

To develop a direct-infusion ESI-MS platform for shotgun lipidomics of human plasma sample, that achieves high-throughput, high-sensitivity, and robustness using Newomics® M3 emitters interfaced with a triple Quadrupole mass spectrometer.

## Introduction

Shotgun lipidomics by ESI-MS is widely used for global analysis of lipids in biological samples such as human plasma. Direct-infusion ESI-MS by nanoflow (flow rate < 1  $\mu\text{L}/\text{min}$ ) is routinely performed to achieve high sensitivity for identifying and quantifying hundreds of lipid species in small volumes of samples. Current single-nozzle emitter technologies including PicoTip<sup>1,2,3,4</sup> and NanoTip<sup>5</sup> have been utilized for shotgun lipidomics at the nanoflow. The Newomics award-winning silicon multinozzle emitters (M3 emitters) split the incoming nanoflow eluent into multiple sub-nanoflows at each nozzle, thereby significantly enhancing the ionization efficiency and reducing the matrix effects for ESI-MS<sup>6,7,8,9,10</sup>. In this Application Note, we demonstrate the optimization of sensitivity, throughput, and robustness for shotgun lipidomics of human plasma samples, using Newomics® M3 emitters interfaced with a Thermo Fisher TSQ Quantiva™ mass spectrometer. The identification and quantification of 577 lipid species from 12 classes in 5  $\mu\text{L}$  human plasma were achieved using precursor and neutral loss modes in a 30-min scan.

## Key Words

M3 emitter, Multinozzle,  
Direct infusion, ESI-MS,  
Nanospray, Shotgun lipidomics,  
Human plasma,  
TSQ Quantiva Triple Quadrupole MS

## Methods

### 1. Sample preparation

Pooled human plasma samples (NIST SRM1950) were thawed at room temperature before lipid extraction. The solvents including chloroform, isopropanol, water, and methanol were HPLC-grade and purchased from Sigma-Aldrich (St. Louis, MO). The lipid standards other than free fatty acids were purchased from Avanti Polar Lipid (Alabaster, AL). The isotope-labeled free fatty acids, DHA d5, EPA d5, and AA d11, were purchased from Sigma-Aldrich (St. Louis, MO). The other lipid internal standards are PC 33:1 d7, PC 48:2, PE 33:1 d7, PE 4ME 16:0 diether, PG 33:1 d7, PG 4ME 16:0 diether, PI 33:1 d7, LPC 18:1 d7, LPC 14:0, LPE 18:1 d7, CE 18:1d7, DAG 33:1 d7, TAG 48:1 d7, SM 36:2 d9, SM 30:1, and Cer 17:0. Total lipids were extracted from 5  $\mu$ L plasma using a modified Bligh & Dyer method<sup>11</sup>. The extracted total lipids were dried under gas nitrogen and reconstituted in acetonitrile: isopropanol: water (65:30:5, v/v/v, 7.5 mM ammonium formate) for direct infusion ESI-MS.

### 2. Mass spectrometry analysis

ESI-MS/MS mass spectra were acquired on a Thermo Fisher TSQ Quantiva triple quadrupole mass spectrometer interfaced with an UltiMate 3000 RSLC nano UPLC system (Thermo Fisher Scientific). The conditions for mass spectrometer settings were listed on **Table 1**. The 12.5  $\mu$ L lipid

sample was continuously infused into the mass spectrometer at a flow rate of 500 nL/min through a New Objective PicoTip™ emitter (single nozzle, 30  $\mu$ m ID, catalog # FS360-50-30-D20), a Thermo Fisher stainless steel emitter (NanoTip, single nozzle, 30  $\mu$ m ID, catalog # NC0379152), or a Newomics® M3 emitter (8 nozzles, 10  $\mu$ m ID, catalog # E8N10MU01). The infusing solvent consisted of acetonitrile: isopropanol: water (65:30:5, v/v/v) and 7.5 mM ammonium formate, which was the same as the sample solvent. The scan rate was 1000 Da/sec. Collision induced dissociation (CID) gas was set at 0.5 mTorr. No source fragmentation voltage was applied. The ion transfer tube temperature was set at 300 °C and the electrospray voltage was optimized and set at 1.8 kV for PicoTip, 2.2 kV for NanoTip, and 3.5 kV for M3 emitter, respectively. Different classes of lipids were scanned at different modes. Precursor ion scan (Prec184) was applied to detect PC, LPC, and SM. Neutral loss scan (NL141) was applied to detect PE and LPE. The PI, PG, CE, and Cer were detected at NL277, NL189, Prec369, and Prec264, respectively. The collision energy applied in each scan mode was 40 V (Prec184), 30 V (NL141), 36 V (NL277), 24 V (NL189), 36 V (Prec369), and 22 V (Prec264), respectively.

**Table 1:** Mass spectrometer parameters for shotgun lipidomics using three types of emitters

	PicoTip	M3 emitter	NanoTip
Flow rate (nL/min)	500		
Mass spectrometer	Thermo TSQ Quantiva		
ID and nozzle number	30 $\mu$ m, 1 nozzle	10 $\mu$ m, 8 nozzles	30 $\mu$ m, 1 nozzle
Scan rate	1000 Da/sec		
Collision induced dissociation gas	0.5 mTorr		
Electrospray voltage	1.8 kV	3.5 kV	2.2 kV
Ion transfer tube temperature	300 °C		
Gas flow (L/min)	0		
Injection mode and volume	User defined, 15 $\mu$ L out of 20 $\mu$ L loop		
Sample solvent	Plasma lipid extract in acetonitrile:Isopropanol:Water (65:30:5, v/v/v) and 7.5 mM ammonium formate		
Run time	30 min		
Infusing solvent or mobile phase	Acetonitrile:Isopropanol:Water (65:30:5, v/v/v) and 7.5 mM ammonium formate		

DAG and TAG were scanned by the neutral loss of fatty acyl residues. These multiple neutral loss scans of fragmentation enabled the quantification of isobaric TAG species<sup>12, 13, 14</sup>. The collision energy for the loss of fatty acyl residues was 32 V. The mass analyzers were adjusted to achieve a resolution of 0.7 atomic mass (am) unit full width at the half height.

### 3. Lipid identification and quantitation

The lipid internal standards and parameters of each scan for shotgun lipidomics are listed on **Table 2**. Raw MS data was processed using Thermo Xcalibur 3.0 and MSFileReader software (Thermo Fisher Scientific). Each lipid class was scanned at a specific mode and one continuum spectrum was acquired. At Prec 184 mode specifically for PC, there were 29 continuum scans. Using the Xcalibur 3.0 software, the 29 continuum scans of Prec 184 were averaged, smoothed, and the peak intensities were determined in the centroid mode at each mass-to-charge ratio (m/z). The data list of mass-to-charge ratios and the corresponding peak intensities were copied

and pasted into a Microsoft Excel sheet. Besides the Prec 184, the 88 continuum scans of NL 141 were also averaged, smoothed, and the peak intensities were pasted into a Microsoft Excel sheet. Similar procedures were performed for NL 277 (79 continuum scans), NL 189 (130 continuum scans), Prec 369 (40 continuum scans), Prec 264 (196 continuum scans), and neutral loss of each fatty residue (29 continuum scans). Free fatty acids, free DHA, free EPA, and free AA, were detected by selected ion monitoring (SRM) and quantified by the corresponding isotope-labeled internal standard. A software tool developed in house at Newomics was used to automatically average and smooth the continuum scans at each scan mode and transfer peak intensities into separate Microsoft Excel sheets. Identification of the peaks of interest and calculation of lipid species amounts were performed using consoles in Microsoft Excel. Corrections for overlap of isotopic variants were applied. The lipids in each class were quantified using the internal standards of that class, respectively. The lipids chosen as the internal standards were not present endogenously in measurable quantities in plasma samples.


**Table 2:** Parameters for each of the 12 lipid classes in our shotgun lipidomics platform

	Internal standard	Scan mode	Collision energy (V)
Phosphatidylcholine (PC)	PC 33:1 d7; PC 48:2	Prec 184, positive	40 V
Phosphatidylethanolamine (PE)	PE 33:1 d7; PE 4ME 16:0 diether	NL 141, positive	30 V
Phosphatidylglycerol (PG)	PG 33:1 d7; PG 4ME 16:0 diether	NL 189, positive	24 V
Phosphatidylinositol (PI)	PI 33:1 d7	NL 277, positive	36 V
Lysophosphatidylcholine (LPC)	LPC 18:1 d7; LPC 14:0	Prec 184, positive	40 V
Lysophosphatidylethanolamine (LPE)	LPE 18:1 d7	NL 141, positive	30 V
Cholesteryl Ester (CE)	CE 18:1 d7	Prec 369, positive	36 V
Diacylglycerol (DAG)	DAG 33:1 d7	Multiple NL, positive	32 V
Triacylglycerol (TAG)	TAG 48:1 d7	Multiple NL, positive	32 V
Sphingomyelin (SM)	SM 36:2 d9; SM 30:1	Prec 184, positive	40 V
Ceramide (Cer)	Cer 17:0	Prec 264, positive	22 V
Free fatty acid (FFA)	DHA d5, EPA d5, AA d11	SRM, positive	12 V

## Results and Discussion

### 1. A shotgun lipidomics platform for human plasma


We established a shotgun lipidomics method by direct infusion ESI-MS/MS analysis of pooled plasma samples. The 12.5  $\mu$ L lipid sample was continuously infused into the mass spectrometer at a flow rate of 500 nL/min (**Figure 1A**). During the 30-min direct infusion, the mass spectra of 12 lipid classes, including charged and neutral lipid classes, were consecutively acquired. The 12 classes of lipids were identified and quantified, including phosphatidylcholine (PC), lysophosphatidylcholine (LPC), phosphatidylethanolamine (PE), lysophosphatidylethanolamine (LPE), phosphatidylinositol (PI), phosphatidylglycerol (PG), sphingomyelin (SM), ceramide (Cer), diacylglycerol (DAG), cholesterol ester (CE), triacylglycerol (TAG), and free fatty acid (FFA). Detailed scan parameters for each lipid class were displayed on **Table 2**. As a representative, neutral loss mode was applied to scan PE and the acquired spectra were shown (**Figure 1B**). Our shotgun lipidomics platform quantified 577 lipid species in human plasma, including 61 PC species, 25 LPC species, 59 PE species, 17 LPE species, 12 PI species, 21 PG species, 41 SM species, 14 Ceramide species, 10 DAG species, 29 CE species, 185 TAG species, and 3 free fatty acid species (**Figure 1C**).




**Figure 1.** Shotgun lipidomics of pooled human plasma showing identification and quantification of 577 lipid species. **A.** TIC of lipid spectra acquired by direct infusion ESI-MS of one sample. Each lipid class was scanned in a separate segment. **B.** Spectra of the neutral loss 141 scan for PE species. **C.** Number of lipid species in each of the 12 lipid classes that were identified and quantified using our shotgun lipidomics platform.

## 2. Sensitivity increase for lipid analysis using M3 emitters

We examined the sensitivity gain of M3 multinozzle emitters over single-nozzle emitters. The identical lipid sample was infused into mass spectrometer at the same flow rate (500 nL/min), for three different types of emitters. The M3 emitter detected significantly higher signal intensity in multiple lipid classes. As an example, **Figure 2** compares signal intensity for three representative PE species (high to low abundance) in 10 consecutive runs using three types of emitters, respectively. On average, M3 emitters achieved 2.5-fold increase in sensitivity compared to PicoTip and NanoTip emitters. Among the 59 PE species, a majority of them displayed over 2-fold higher signal using the M3 emitters than those using PicoTip and NanoTip emitters, respectively. **Figure 3** shows the sensitivity gain of the M3 emitter over PicoTip and NanoTip for the top 25 PE species in human plasma.




**Figure 2.** Representative TIC signals of PE species obtained using different type of emitters. The exactly same plasma lipid samples were infused into the mass spectrometer and the signals were acquired. We performed 10 consecutive injections using M3 emitter, NanoTip, and PicoTip, respectively. **A.** The most abundant PE species detected in plasma. **B.** The medium abundant PE species in plasma. **C.** The least abundant species in plasma.



**Figure 3.** Sensitivity gain of ESI-MS using M3 emitter compared to single-nozzle emitters on the most abundant 25 PE species in human plasma. The signals of PicoTip and NanoTip emitters were used as denominators. On average, 2-fold increase in sensitivity was observed for M3 emitters.

### 3. Improvement in reproducibility and LOD for lipid analysis using M3 emitters

We compared the CV for 10 consecutive runs for the most abundant 10 PE species using three types of emitters (**Figure 4A**). On average, the M3 emitter achieved ~5% CV, comparable to those using PicoTip, while significantly outperforming those using NanoTip (~7.5% CV). We evaluated the LOD using PE 33:1 d7, an isotope-labelled internal standard. The M3 emitter displayed the lowest LOD (S/N=3) of 0.17 pmole/ $\mu$ L (**Figure 4B**), significantly outperforming both PicoTip and NanoTip. This is consistent with the higher sensitivity gained by M3 emitters (**Figure 2 and 3**).



**Figure 4.** Reproducibility and detection limit of PE species by ESI-MS using M3, NanoTip, and PicoTip emitters. **A.** CV of the signals of the 10 most abundant PE species (n=10). **B.** LOD for one representative PE species PE 33:1 d7 for three types of emitters.

## Conclusions

We have established a new ESI-MS platform for shotgun lipidomics studies of human plasma samples. By interfacing to a triple Quadrupole mass spectrometer, we demonstrated that Newomics® M3 emitters could achieve high throughput, high sensitivity, and robustness for shotgun lipidomics. For most of the lipid classes in human plasma, our M3 emitters significantly outperform the conventional single-nozzle emitters in sensitivity, robustness, and LOD for their identification and quantitation using ESI-MS/MS analysis.

## References

1. Horn PJ, Ledbetter NN, James CN, Hoffman WD, Case CR, Verbeck GF, Chapman KD, (2010), *Journal of Biological Chemistry*, 286, 3298-3302.
2. Holzer M, Birner-Gruenberger R, Stojakovic T, El-Gamal D, Binder V, Wadsack C, Heinemann A, Marsche G, (2011), *Journal of the American Society of Nephrology*, 22, 1631-1641.
3. Phelps MS, Verbeck GF, (2020), *Single Cell Metabolism*, 19-30.
4. Sivedoroudi P, Bennike TB, Kastaniegaard K, Talebpour M, Ghassempour A, Stensballe A, (2019), *Journal of Proteomics*, 203, 103373.
5. Schilcher I, Ledinski G, Radulovic S, Hallström S, Eichmann T, Madl T, Zhang F, Leitinger G, Kolb-Lenz D, Darnhofer B, Birner-Gruenberger R, Wadsack C, Kratky D, Marsche G, Frank S, Cvirn G, (2019), *Biochim Biophys Acta - Molecular and Cell Biology of Lipids*, 1864(10), 1363-1374.
6. Kim W., Guo M., Yang P., Wang D., (2007), *Analytical Chemistry*, 79 (10), 3703-3707.
7. Mao P., Wang H.T., Yang P., Wang D., (2011), *Analytical Chemistry*, 83(15), 6282-6289.
8. Mao P., Gomez-Sjoberg R., Wang D., (2013), *Analytical Chemistry*, 85(2), 816-819.
9. Newomics Application Note 1:  
<https://www.newomics.com/wp-content/uploads/2019/05/Newomics-Application-Note-1.pdf>
10. Newomics Application Note 2:  
<https://www.newomics.com/wp-content/uploads/2019/05/Newomics-Application-Note-2.pdf>
11. Bligh and Dyer, 1959, A rapid method of total lipid extraction and purification.
12. Han, X.; Gross, R. W. (2001), *Anal Biochem* 2001, 295, 88-100.
13. Li, M.; Butka, E.; Wang, X. (2014), *Sci Rep* 2014, 4, 6581.
14. Li, M.; Baughman, E.; Roth, M. R.; Han, X.; Welti, R.; Wang, X. (2014), *Plant J* 2014, 77, 160-172.

## Ordering Information

Catalog #	Product
SKIT-T01	Newomics® Starter Kit for Thermo Fisher Nanospray Flex Ion Sources
E8N10MU01	Newomics® M3 Emitters, 10 µm I.D. – 8 nozzles

Newomics Inc.  
804 Heinz Ave, STE 150, Berkeley, CA 94710  
To learn more, please contact: [sales@newomics.com](mailto:sales@newomics.com)  
[www.newomics.com](http://www.newomics.com)

© 2019 Newomics Inc. All rights reserved. The information is presented as an example of the capabilities of Newomics products. Specifications, terms, and pricing are subject to change.

